

Laparoscopic Port-Site Fascial Closure Device Review and Catalog

Prindeze N, Foshag A. Chiba S, Pechman D, Renton D

Review of evidence for port site fascial closure

A number of factors have been identified which are significantly correlated with hernia development and reoccurrence. These include non-modifiable factors such as age and male gender, and modifiable ones such as obesity, diabetes, immunosuppression, smoking and COPD^{1,2}. It is important patients and providers still address these underlying risk factors even though just the use of a laparoscopic port site is not traditionally considered a hernia.

Laparoscopic surgery has evolved significantly since its introduction nearly forty years ago. Understanding nuances which contribute to port site hernia occurrence can therefore be challenging as technique and technology have considerably changed. Factors, however which have been consistently associated with port-site hernia development include trochar size, abdominal location and port-site enlargement.

Port size

Perhaps one of the most studied factors in port site hernias is trochar size. A 1997 prospective study identified 5mm trochars to be associated with only one of 16 hernias (6.3%) in 870 patients³. Similarly, a 2010 systematic review identified only 2 cases (2.0%) from 5mm ports of 99 hernias⁴, and this trend continues through more recently published work as well; 2.7% in a 2011 systematic review⁵, and 3% in a 2017 retrospective review⁶. All other hernias in these studies occurred with ports 10mm or greater. A 2020 systematic review by Gutierrez compared port-site closure vs. none for both 5mm and 10mm trochars and identified a herniation rate of 0% and 0.08% for 5mm ports and 0.14% and 0.2% for 10mm ports⁷.

Unequivocally single port systems have been associated with significantly greater rates of port-site hernias with a 2021 single center trial reporting an 8% incidence at 2 years⁸, and a 2018 review of 13 studies reporting a 5.2% incidence at a median follow-up of 4 months⁹.

Abdominal location

Several factors regarding trochar location on the abdominal wall have been associated with increased incidence of port-site hernia. Overall hernias occur most commonly at the umbilical port site but have been observed in all port-site locations.

A 2010 review by Bunting et.al. identified 89% of 99 hernias to be at the umbilicus with the remaining 11% upper abdominal⁴. A 2017 review by Lambertz et. al. identified a more even distribution with 44% of 54 hernias umbilical or peri-umbilical, 35% epigastric and 20% lateral abdominal wall⁶. A 2020 review by Gutierrez calculated incidence by site with the highest of 1.2% at the umbilicus, extraumbilical at 0.33% and subxiphoid 0.60%⁷.

One repeatedly identified factor associated with port-site hernia is midline location. In a 2011 systematic review Owens et. al. observed that 57% of port-site hernias occurred in the midline⁵

while Gutierrez et. al. determined an overall midline port-site hernia incidence of 0.71% as compared to off-midline at 0.28%, a 2.5 fold difference⁷.

Port-site extension

Port-site extension or enlargement for specimen extraction was one of the earlier identified factors associated with hernia development. Nassar in 1997 found extension to have been involved in 12 of 16 (75%) port-site hernias³. Ki et. al in 2019 identified extension to have occurred in 10 of 18 (56%) hernias¹⁰, and Bunting et. al in 2010 identified extension to have been a significant risk factor in a majority of included studies⁴.

Other factors

In their 1997 prospective study, Nassar et. al. identified pre-existing umbilical hernia as a significant risk factor for post-operative port-site hernia³. Of 870 laparoscopic cholecystectomies, 12% presented with an incidental finding of umbilical hernia which was subsequently used as a port-site and closed primarily with suture. However, of sixteen port-site hernias identified in their study, four (25%) had hernias pre-operatively. Bunting et. al similarly identified pre-existing umbilical hernia as a significant risk factor for port-site hernia development in their 2010 systematic review⁴.

Other risk factors suggested for port-site hernia include: Older age^{4,10}, higher BMI¹⁰, sarcopenia¹⁰, length of surgery¹⁰, presence of a rectus diastasis¹⁰, port-site infection³⁻⁵ and male gender^{3,4}.

Devices

Product	Entry	Disposable	Unprotected sharp	Requires a trocar	Fits through a trocar	Functional length	Width
Carter-Thomason CloseSure®	1998	Yes	Yes	No	No	6 – 10 cm	**
Carter-Thomason II™	2012	Yes	Yes	No	Yes	11 – 15 cm	**
Efx Classic™	2010	Yes	Yes	No	No	5.7 - 10.1 cm	1.3 – 1.7cm
Efx Shield®	2015	Yes	No	No	No	6.1 cm	1.3 cm
CrossBow™ and Longbow™	2015	Yes	Yes	No	No	7cm 9cm	1.4cm
EZ Stitch™	2020	Yes	No	No	No	17cm	1.4cm
Lapro-Shark™	2016	Yes	No	No	No	7.6 cm	1.2 - 2.9 cm
SutureClose™	2016	Yes	Yes	No	No	6 cm	1.4 cm
VersaOne™	2017	Partially	Yes	Yes	Yes (Limited)	5.1cm	1.2cm
NeoClose®	2018	Yes	Yes	No	No	**	**


M-close™	2019	Yes	No	No	Yes	10.2 cm	1 cm
----------	------	-----	----	----	-----	---------	------

**Manufacturer prefers not to disclose/Proprietary

Carter-Thomason CloseSure System®

Manufacturer: Cooper Surgical®

Entry Date: 1998



Images reprinted with permission of Cooper Surgical®

Mechanism: Suture passer needle paired with cone tip guide

Components:

1. Needle guide
2. Suture retriever needle

Disposability: Disposable

Unprotected Sharp: Yes

Requires a trochar: No

Fits through a trochar: No

Fascial bite: Variable based on depth

Functional dimensions:

Length (Bottom of the flange to the needle exit sites)

5mm Guide:

10/12 Guide: 6 cm

15mm Guide: 7 cm

XL 10/12 Guide: 10 cm

XL 15mm Guide: 10 cm

Width (Maximal cone diameter): Manufacturer prefers not to disclose


Cooper Surgical CloseSure System ([Documentation Link](#))

Cooper Surgical Product Comparison ([Documentation Link](#))

Carter-Thomason® II Port Closure System

Manufacturer: Cooper Surgical®

Entry Date: 2012



Images reprinted with permission of Cooper Surgical®

Mechanism: Suture passer needle paired with stick guide with distal collapsable non-needle-obstructing wings

Components:

1. Needle guide
2. Suture retriever needle

Disposability: Disposable

Unprotected Sharp: Yes

Requires a trochar: No

Fits through a trochar: Yes

Fascial bite: Variable based on depth

Functional dimensions:

Length (Top of the anchoring wings to needle entry point)

10/12 System: 11 cm

15mm System: 11 cm

XL 10/12 System: 15 cm

XL 15mm System: 14 cm

Width (Maximal width of the shaft): Manufacturer prefers not to disclose

Cooper Surgical Carter Thomason II ([Documentation Link](#))

Cooper Surgical Product Comparison ([Documentation Link](#))

EFx Classic™ Fascial Closure System

Manufacturer: Teleflex™ Weck®

Entry Date: 2010



Images reprinted with permission of Teleflex™ Weck®

Mechanism: Suture passer needle paired with cone tip guide

Components:

1. Needle guide
2. Suture retriever needle

Disposability: Disposable with optional reusable Berci needle passer

Unprotected Sharp: Yes

Requires a trochar: No

Fits through a trochar: No

Fascial bite: Variable based on depth

Functional dimensions:

Length (Bottom of the flange to the needle exit sites):

10/12 Guide: 5.7 cm

15mm Guide: 10.1 cm

Width (Maximal cone diameter):

10/12 Guide: 1.3 cm


15mm Guide: 1.7 cm

Teleflex Weck EFx Shield and EFx Classic ([Documentation Link](#))

EFx Shield® Fascial Closure System

Manufacturer: Teleflex™ Weck®

Entry Date: 2015



Images reprinted with permission of Teleflex™ Weck®

Mechanism: Suture passer needle paired with stick guide with distal collapsable needle-obstructing wings

Components:


1. Needle guide
2. Suture retriever needle

Disposability: Disposable**Unprotected Sharp:** No**Requires a trochar:** No**Fits through a trochar:** No**Fascial bite:** 1cm (Guide shaft to needle entry site on wing)**Functional dimensions:**

Length: 6.1 cm (Top of wings to tapered neck below suture entry)

Width: 1.3 cm (Shaft diameter)

Teleflex Weck EFX Shield and EFX Classic ([Documentation Link](#))**CrossBow™ and LongBow™ Closure Device****Manufacturer:** Suture Ease**Entry Date:** 2015



Images reprinted with permission of Suture Ease

Mechanism: Suture passer needle paired with cone tip guide with retractable snare loops

Components:

1. Needle guide
2. Suture retriever needle

Disposability: Disposable

Unprotected Sharp: Yes

Requires a trochar: No

Fits through a trochar: No

Fascial bite: Variable based on depth (Indicated with marks along shaft)

Functional dimensions:

Crossbow

Length: 7cm (Bottom of the flange to the snares)

Width: 1.4cm (Maximal width of the cone)

Longbow

Length: 9cm (Bottom of the flange to the snares)

Width: 1.4cm (Maximal width of the cone)


Suture Ease CrossBow ([Documentation Link](#))

Suture Ease CrossBow Bariatric ([Documentation Link](#))

EZ Stitch™ Port Site Closure System

Manufacturer: Suture Ease

Entry Date: 2020



Images reprinted with permission of Suture Ease

Mechanism: Stick-style needle guide with "C"-shaped guide slot and suture passer

Components:

1. Needle guide
2. Suture retriever needle

Disposability: Disposable

Unprotected Sharp: No

Requires a trocar: No

Fits through a trocar: No

Fascial bite: 7.5mm (Needle to back of "C" slot)

Functional dimensions:

Length: 17cm (skin-contacting flange to the bottom of the "C"-shaped fascial catch)



Width: 1.4cm diameter

Suture Ease EZ Stitch ([Documentation Link](#))

Lapro-Shark™ Laparoscopic Port Closure Device

Manufacturer: Brainchild Surgical Devices

Entry Date: 2016



Images reprinted with permission of Brainchild Surgical Devices

Mechanism: "C"-shaped guide slot with needle guide

Components:

1. Needle guide
2. Suture retriever needle

Disposability: Disposable

Unprotected Sharp: No

Requires a trocar: No

Fits through a trocar: No

Fascial bite: 1cm (Needle to back of “C” slot)

Functional dimensions:

Length: 7.6cm (Bottom of the "C" to the entry point of the needle)



Width: 1.2 - 2.9cm (Tapered)

Brainchild Surgical LaproShark (Documentation Link)

SutureClose™ Laparoscopic Port Closure System

Manufacturer: DeRoyal®

Entry Date: 2016



Images reprinted with permission of DeRoyal®

Mechanism: Suture passer needle paired with cone tip guide

Components:

1. Needle guide
2. Suture retriever needle

Disposability: Disposable

Unprotected Sharp: Yes

Requires a trochar: No

Fits through a trochar: No

Fascial bite: Variable based on depth

Functional dimensions:

Length: 6 cm (Bottom of the flange to the needle exit sites)

Width: 1.4 cm (Maximal cone diameter)


DeRoyal SutureClose ([Documentation Link](#))

VersaOne™ Fascial Closure System

Manufacturer: Medtronic Covidien

Entry Date: 2017





Images reprinted with permission of Medtronic Covidien

Mechanism: Suture passer needle guide obturator placed through existing trochar

Components:

1. Trochar
2. Obturator (Bladed or Optical)
3. Guide
4. Suture passer (sold separately)

Disposability: Disposable, however requires the reusable suture passer

Unprotected Sharp: Yes

Requires a trochar: Yes (Can be used as a needle guide without the trochar)

Fits through a trochar: Yes (Only designated VersaOne trochar)

Fascial bite: Variable based on depth

Functional dimensions:

Length: 5.1 cm (Bottom of the flange to the needle exit sites)




Width: 1.2 cm (Trochar diameter)

VersaOne Fascial Closure System ([Documentation Link](#))

NeoClose® AnchorGuard

Manufacturer: Sharp Fluidics®

Entry Date: 2018



Images reprinted with permission of Sharp Fluidics®

Mechanism: Suture passer needle with pre-loaded anchors and suture paired with cone tip guide

Components:

1. Needle guide
2. Suture passers with pre-loaded PDS anchors and sutures

Disposability: Disposable

Unprotected Sharp: Yes

Requires a trochar: No

Fits through a trochar: No



Fascial bite: Variable based on depth

Functional dimensions:

Length: 9.9 cm (Total length)
 Width (Maximal cone diameter): Manufacturer prefers not to disclose – Recommended ranges: Small guide (5-12mm), Large guide (8-15mm)
 Sharp Fluidics NeoClose (Documentation Link)

M-Close™ Kit

Manufacturer: New Wave Endo-Surgery™
Entry Date: 2019



Images reprinted with permission of New Wave Endo-Surgery™

Mechanism: Stick device with deployable anvil and needles that a wire and suture are passed through. Local anesthetic can also be used through the needles

Components:

1. Stick with anvil and needles
2. Wire
3. Wire passer
4. Kit also includes separate nerve block needle

Disposability: Disposable

Unprotected Sharp: No

Requires a trocar: No

Fits through a trocar: Yes

For ports <12mm, the trocar is removed and the M-close is inserted through the skin
 Will fit through bariatric trochars up to 150mm in length

Fascial bite: 1cm (Shaft to needle entry site on T-bar)

Functional dimensions: 1cm

Length: 4 inches (10.2cm) (Top of the T-bar to the slider button)

Width: 1.2 cm

New Wave Endo M-close ([Documentation Link](#))

References

- (1) Sørensen, L. T.; Hemmingsen, U. B.; Kirkeby, L. T.; Kallehave, F.; Jørgensen, L. N. Smoking Is a Risk Factor for Incisional Hernia. *Arch. Surg.* **2005**, *140* (2), 119–123. <https://doi.org/10.1001/archsurg.140.2.119>.
- (2) Bhardwaj, P.; Huayllani, M. T.; Olson, M. A.; Janis, J. E. Year-Over-Year Ventral Hernia Recurrence Rates and Risk Factors. *JAMA Surg.* **2024**, *159* (6), 651–658. <https://doi.org/10.1001/jamasurg.2024.0233>.
- (3) Nassar, A. H.; Ashkar, K. A.; Rashed, A. A.; Abdulmoneum, M. G. Laparoscopic Cholecystectomy and the Umbilicus. *Br. J. Surg.* **1997**, *84* (5), 630–633.
- (4) Bunting, D. M. Port-Site Hernia Following Laparoscopic Cholecystectomy. *JSLs* **2010**, *14* (4), 490–497. <https://doi.org/10.4293/108680810X12924466007728>.
- (5) Owens, M.; Barry, M.; Janjua, A. Z.; Winter, D. C. A Systematic Review of Laparoscopic Port Site Hernias in Gastrointestinal Surgery. *Surg. J. R. Coll. Surg. Edinb. Irel.* **2011**, *9* (4), 218–224. <https://doi.org/10.1016/j.surge.2011.01.003>.
- (6) Lambertz, A.; Stüben, B. O.; Bock, B.; Eickhoff, R.; Kroh, A.; Klink, C. D.; Neumann, U. P.; Krones, C. J. Port-Site Incisional Hernia - A Case Series of 54 Patients. *Ann. Med. Surg.* **2017**, *14*, 8–11. <https://doi.org/10.1016/j.amsu.2017.01.001>.
- (7) Gutierrez, M.; Stuparich, M.; Behbehani, S.; Nahas, S. Does Closure of Fascia, Type, and Location of Trocar Influence Occurrence of Port Site Hernias? A Literature Review. *Surg. Endosc.* **2020**, *34* (12), 5250–5258. <https://doi.org/10.1007/s00464-020-07826-8>.
- (8) Valcarenghi, J.; Hernigou, J.; Apicella, G.; Clegg, E.; Rousie, M.; Chasse, E. Long-Term Follow-up of the Incisional Hernia Rate after Single-Incision Laparoscopic Cholecystectomy: A Prospective Observational Study. *Acta Chir. Belg.* **2021**, *121* (5), 320–326. <https://doi.org/10.1080/00015458.2020.1765673>.
- (9) Migliore, M.; Arezzo, A.; Arolfo, S.; Passera, R.; Morino, M. Safety of Single-Incision Robotic Cholecystectomy for Benign Gallbladder Disease: A Systematic Review. *Surg. Endosc.* **2018**, *32* (12), 4716–4727. <https://doi.org/10.1007/s00464-018-6300-3>.
- (10) Ki, H. J.; Park, J. B.; Sul, J.-Y. Umbilical Port Site Hernia and Diastasis Recti. *J. Minim. Invasive Surg.* **2020**, *23* (2), 80–85. <https://doi.org/10.7602/jmis.2020.23.2.80>.